



St Andrew's First Aid

# Defibrillator checklist

We've put together a checklist to help you ensure your defibrillator is always ready for use.

Defibrillators self-test on a daily, weekly or monthly basis – so your device will signal, typically with a flashing light or audible alert, if there is a problem, for example, pads not correctly connected, or low battery.



- 1 Set up a regular inspection procedure for your defibrillator**
- 2 Record your inspections below**
- 3 Record the expiration date of pads and batteries**
- 4 Reorder batteries and pads before they expire**
- 5 Follow the guidance in your product manual when cleaning your defibrillator**

Pads expiration date

Batteries expiration date

| Defibrillator ready to use? | Date | Initial |
|-----------------------------|------|---------|
|                             |      |         |
|                             |      |         |
|                             |      |         |
|                             |      |         |
|                             |      |         |
|                             |      |         |
|                             |      |         |
|                             |      |         |
|                             |      |         |

Need new pads or batteries? Order now at [www.firstaid.org.uk/supplies](http://www.firstaid.org.uk/supplies) or call 0300 4 666 999 option 1/option 3 (we are available Monday – Friday 9am to 5pm)